

# Thrombin-like enzyme crotalase Antibody



PACO35090

---

## Product Information

**Size:**

50ug

**Reactivity:**

Crotalus adamanteus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Thrombin-like snake venom protein that release fibrinopeptide A from fibrinogen (FGA). Shows both kinin-releasing and coagulant activities.

**Gene ID:**

**Uniprot**

F8S114

**Synonyms:**

Thrombin-like enzyme crotalase (SVTLE) (EC 3.4.21.74) (Fibrinogen-clotting enzyme) (Snake venom serine protease 2) (SVSP)

**Immunogen:**

Recombinant Crotalus adamanteus Thrombin-like enzyme crotalase protein (25-262AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A