

PACO35118

Product Information

Size:

50ug

Reactivity:

Bacillus subtilis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Hydrolysis of terminal, non-reducing alpha-D-galactose residues in alpha-D-galactosides, including galactose oligosaccharides, galactomannans and galactolipids.

Gene ID:

meIA

Uniprot

O34645

Synonyms:

Alpha-galactosidase (EC 3.2.1.22) (Melibiase), meIA

Immunogen:

Recombinant Bacillus subtilis Alpha-galactosidase protein (1-432AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A