

# X Antibody



PACO35134

---

## Product Information

**Size:**

50ug

**Reactivity:**

Hepatitis B virus genotype D

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

X

**Uniprot**

O93195

**Synonyms:**

Protein X (HBx) (Peptide X) (pX), X

**Immunogen:**

Recombinant Hepatitis B virus genotype D Protein X protein (1-154AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A