

CSNK2A2 Antibody, Biotin conjugated



PACO35149

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	CSNK2A2
Source:	Uniprot
Rabbit	P19784
Isotype:	Synonyms:
IgG	Casein kinase II subunit alpha' (CK II alpha') (EC 2.7.11.1), CSNK2A2, CK2A2
Applications:	Immunogen:
ELISA	Recombinant Human Casein kinase II subunit alpha' protein (8-350AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A