

# dsbA Antibody



PACO35162

---

## Product Information

**Size:**

50ug

**Reactivity:**

Pseudomonas syringae pv. tomato

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in disulfide-bond formation. Acts by transferring its disulfide bond to other proteins.

**Gene ID:**

dsbA

**Uniprot**

O52376

**Synonyms:**

Thiol: disulfide interchange protein DsbA, dsbA

**Immunogen:**

Recombinant Pseudomonas syringae pv. tomato Thiol: disulfide interchange protein DsbA protein (23-214AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A