

PACO35174

Product Information

Size:

50ug

Reactivity:

Bacillus thuringiensis subsp. morrisoni

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes colloidosmotic lysis by binding to the midgut epithelial cells of insects.

Gene ID:

cry1Fb

Uniprot

O66377

Synonyms:

Pesticidal crystal protein Cry1Fb (132 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIF(b)), cry1Fb, cryIF(b) cryINA67-1

Immunogen:

Recombinant Bacillus thuringiensis subsp. morrisoni Pesticidal crystal protein Cry1Fb protein (984-1169AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A