

rpiA Antibody, HRP conjugated



PACO35195

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the reversible conversion of ribose-5-phosphate to ribulose 5-phosphate.

Gene ID:

rpiA

Uniprot

P0A7Z0

Synonyms:

Ribose-5-phosphate isomerase A (EC 5.3.1.6) (Phosphoriboisomerase A) (PRI), rpiA, ygfC

Immunogen:

Recombinant Escherichia coli Ribose-5-phosphate isomerase A protein (1-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A