

C2orf80 Antibody



PACO35222

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C2orf80

Uniprot

Q0P641

Synonyms:

Uncharacterized protein C2orf80, C2orf80

Immunogen:

Recombinant Human Uncharacterized protein C2orf80 protein (1-193AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A