

# lepB Antibody



PACO35226

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mycobacterium tuberculosis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

lepB

**Uniprot**

P9WKA1

**Synonyms:**

Signal peptidase I (SPase I) (EC 3.4.21.89) (Leader peptidase I), lepB

**Immunogen:**

Recombinant Mycobacterium tuberculosis Signal peptidase I protein (88-294AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A