

ZNF252P-AS1 Antibody



PACO35254

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

ZNF252P-AS1

Uniprot

Q0IIN9

Synonyms:

Putative uncharacterized protein ZNF252P-AS1 (ZNF252P antisense RNA 1) (ZNF252P antisense gene protein 1), ZNF252P-AS1, C8orf77


Immunogen:

Recombinant Human Putative uncharacterized protein ZNF252P-AS1 protein (1-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO35254 at dilution of 1:100.