

ZNF252P-AS1 Antibody, HRP conjugated



PACO35255

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

ZNF252P-AS1

Uniprot

Q0IIN9

Synonyms:

Putative uncharacterized protein ZNF252P-AS1 (ZNF252P antisense RNA 1) (ZNF252P antisense gene protein 1), ZNF252P-AS1, C8orf77

Immunogen:

Recombinant Human Putative uncharacterized protein ZNF252P-AS1 protein (1-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A