

ODF1 Antibody



PACO35310

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the outer dense fibers (ODF) of spermatozoa. ODF are filamentous structures located on the outside of the axoneme in the midpiece and principal piece of the mammalian sperm tail and may help to maintain the passive elastic structures and elastic recoil of the sperm tail.

Gene ID:

ODF1

Uniprot

Q14990

Synonyms:

Outer dense fiber protein 1 (Heat shock protein beta-10) (HspB10), ODF1, HSPB10
ODFP

Immunogen:

Recombinant Human Outer dense fiber protein 1 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A