

SSX2 Antibody



PACO35330

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Could act as a modulator of transcription.

Gene ID:

SSX2

Uniprot

Q16385

Synonyms:

Protein SSX2 (Cancer/testis antigen 5.2) (CT5.2) (Synovial sarcoma, X breakpoint 2) (Tumor antigen HOM-MEL-40), SSX2; SSX2B, SSX2A;

Immunogen:

Recombinant Human Protein SSX2 protein (1-188AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A