

PACO35334

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:**Gene ID:**

NKG7

Uniprot

Q16617

Synonyms:

Protein NKG7 (G-CSF-induced gene 1 protein) (GIG-1 protein) (Granule membrane protein of 17 kDa) (GMP-17) (Natural killer cell protein 7) (p15-TIA-1), NKG7, GIG1


Immunogen:

Recombinant Human Protein NKG7 protein (30-60AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO35334 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).