

# NKG7 Antibody, FITC conjugated



PACO35336

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NKG7

**Uniprot**

Q16617

**Synonyms:**

Protein NKG7 (G-CSF-induced gene 1 protein) (GIG-1 protein) (Granule membrane protein of 17 kDa) (GMP-17) (Natural killer cell protein 7) (p15-TIA-1), NKG7, GIG1

**Immunogen:**

Recombinant Human Protein NKG7 protein (30-60AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A