

C21orf62-AS1 Antibody, HRP conjugated



PACO35355

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C21orf62-AS1

Uniprot

Q17RA5

Synonyms:

Putative uncharacterized protein C21orf62-AS1 (C21orf62 antisense RNA 1) (C21orf62 antisense gene protein 1), C21orf62-AS1, C21orf49

Immunogen:

Recombinant Human Putative uncharacterized protein C21orf62-AS1 protein (1-79AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A