

POLA2 Antibody, Biotin conjugated



PACO35381

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play an essential role at the early stage of chromosomal DNA replication by coupling the polymerase alpha/primase complex to the cellular replication machinery.

Gene ID:

POLA2

Uniprot

Q14181

Synonyms:

DNA polymerase alpha subunit B (DNA polymerase alpha 70 kDa subunit), POLA2

Immunogen:

Recombinant Human DNA polymerase α subunit B protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A