

ZMYM3 Antibody



PACO35382

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:20-1:200

Protein Background:

Plays a role in the regulation of cell morphology and cytoskeletal organization.

Gene ID:

ZMYM3

Uniprot

Q14202

Synonyms:

Zinc finger MYM-type protein 3 (Zinc finger protein 261), ZMYM3, DXS6673E KIAA0385 ZNF261


Immunogen:

Recombinant Human Zinc finger MYM-type protein 3 protein (1-495AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of A549 cells using PACO35382 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).