

# PSG3 Antibody, HRP conjugated



PACO35479

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

PSG3

**Uniprot**

Q16557

**Synonyms:**

Pregnancy-specific beta-1-glycoprotein 3 (PS-beta-G-3) (PSBG-3) (Pregnancy-specific glycoprotein 3) (Carcinoembryonic antigen SG5), PSG3

**Immunogen:**

Recombinant Human Pregnancy-specific beta-1-glycoprotein 3 protein (35-310AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A