

FGL2 Antibody



PACO35586

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play a role in physiologic lymphocyte functions at mucosal sites.

Gene ID:

FGL2

Uniprot

Q14314

Synonyms:

Fibroleukin (Fibrinogen-like protein 2) (pT49), FGL2


Immunogen:

Recombinant Human Fibroleukin protein (24-439AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO35586 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO35586 at dilution of 1:100.