

KRT33B Antibody



PACO35638

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

There are two types of hair/microfibrillar keratin, I (acid, c) and II (neutral to basic).

Gene ID:

KRT33B

Uniprot

Q14525

Synonyms:

Keratin, type I cuticular Ha3-II (Hair keratin, type I Ha3-II) (Keratin-33B) (K33B), KRT33B, HHA3-II HKA3B KRTHA3B


Immunogen:

Recombinant Human Keratin, type I cuticular Ha3-II protein (1-404AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: KRT33B antibody at 1 μ g/ml + MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 47 kDa. Observed band size: 47 kDa.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO35638 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO35638 at dilution of 1:100.