

GABRA6 Antibody, HRP conjugated



PACO35723

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

Gene ID:

GABRA6

Uniprot

Q16445

Synonyms:

Gamma-aminobutyric acid, receptor subunit alpha-6 (GABA(A) receptor subunit alpha-6), GABRA6

Immunogen:

Recombinant Human Gamma-aminobutyric acid, receptor subunit alpha-6 protein (20-242AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A