

KIR2DS1 Antibody



PACO35778

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor on natural killer (NK) cells for HLA-C alleles. Does not inhibit the activity of NK cells.

Gene ID:

KIR2DS1

Uniprot

Q14954

Synonyms:

Killer cell immunoglobulin-like receptor 2DS1 (CD158 antigen-like family member H) (MHC class I NK cell receptor Eb6 Act1) (CD antigen CD158h), KIR2DS1, CD158H

Immunogen:

Recombinant Human Killer cell immunoglobulin-like receptor 2DS1 protein (22-245AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A