

**L-lactate dehydrogenase Antibody, HRP  
conjugated  
PACO35823**



---

**Product Information**

**Size:**

50ug

**Protein Background:**

**Reactivity:**

Plasmodium falciparum

**Gene ID:**

**Source:**

Rabbit

**Uniprot**

Q27743

**Isotype:**

IgG

**Synonyms:**

L-lactate dehydrogenase (EC 1.1.1.27) (LDH-P)

**Applications:**

ELISA

**Immunogen:**

Recombinant Plasmodium falciparum L-lactate dehydrogenase protein (238-316AA).

**Recommended dilutions:**

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A