

L-lactate dehydrogenase Antibody, Biotin conjugated
PACO35825



Product Information

Size:

50ug

Reactivity:

Plasmodium falciparum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

Q27743

Synonyms:

L-lactate dehydrogenase (EC 1.1.1.27) (LDH-P)

Immunogen:

Recombinant Plasmodium falciparum L-lactate dehydrogenase protein (238-316AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A