

IZUMO4 Antibody, FITC conjugated



PACO35832

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IZUMO4

Uniprot

Q1ZYL8

Synonyms:

Izumo sperm-egg fusion protein 4 (Sperm 22 kDa protein c113), IZUMO4, C19orf36

Immunogen:

Recombinant Human Izumo sperm-egg fusion protein 4 protein (16-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A