

PACO35834

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Alpha-toxin binds to the membrane of eukaryotic cells resulting in the release of low-molecular weight molecules and leading to an eventual osmotic lysis. Heptamer oligomerization and pore formation is required for lytic activity.

Gene ID:

hly

Uniprot

Q2G1X0

Synonyms:

Alpha-hemolysin (Alpha-HL) (Alpha-toxin), hly, hla

Immunogen:

Recombinant Staphylococcus aureus Alpha-hemolysin protein (27-319AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A