

# Slc39a13 Antibody, HRP conjugated



PACO35855

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Acts as a zinc-influx transporter.

**Gene ID:**

Slc39a13

**Uniprot**

Q2M1K6

**Synonyms:**

Zinc transporter ZIP13 (Solute carrier family 39 member 13) (Zrt- and Irt-like protein 13) (ZIP-13), Slc39a13, Zip13

**Immunogen:**

Recombinant Rat Zinc transporter ZIP13 protein (130-233AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A