

ZXDC Antibody, Biotin conjugated



PACO35865

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cooperates with CIITA to promote transcription of MHC class I and MHC class II genes.

Gene ID:

ZXDC

Uniprot

Q2QGD7

Synonyms:

Zinc finger protein ZXDC (ZXD-like zinc finger protein), ZXDC, ZXDL

Immunogen:

Recombinant Human Zinc finger protein ZXDC protein (139-416AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A