

# C15orf32 Antibody



PACO35902

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C15orf32

**Uniprot**

Q32M92

**Synonyms:**

Uncharacterized protein C15orf32, C15orf32

**Immunogen:**

Recombinant Human Uncharacterized protein C15orf32 protein (1-165AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A