

C15orf32 Antibody, Biotin conjugated



PACO35905

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

C15orf32

Uniprot

Q32M92

Synonyms:

Uncharacterized protein C15orf32, C15orf32

Immunogen:

Recombinant Human Uncharacterized protein C15orf32 protein (1-165AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A