

ANKRD55 Antibody



PACO35926

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

ANKRD55

Uniprot

Q3KP44

Synonyms:

Ankyrin repeat domain-containing protein 55, ANKRD55


Immunogen:

Recombinant Human Ankyrin repeat domain-containing protein 55 protein (151-326AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO35926 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO35926 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).