

# SLC25A52 Antibody



PACO35958

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC25A52

**Uniprot**

Q3SY17

**Synonyms:**

Solute carrier family 25 member 52 (Mitochondrial carrier triple repeat protein 2), SLC25A52, MCART2

**Immunogen:**

Recombinant Human Solute carrier family 25 member 52 protein (1-33AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A