

GFOD2 Antibody



PACO35966

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200

Protein Background:

Promotes matrix assembly.

Gene ID:

GFOD2

Uniprot

Q3B7J2

Synonyms:

Glucose-fructose oxidoreductase domain-containing protein 2 (EC 1. -. -. -), GFOD2


Immunogen:

Recombinant Human Glucose-fructose oxidoreductase domain-containing protein 2 protein (1-280AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: GFOD2 antibody at 2 μ g/ml. Lane 1: K562 whole cell lysate. Lane 2: HeLa whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 43, 31, 13 kDa. Observed band size: 43 kDa.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO35966 at dilution of 1:100.