

PACO35982

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Probable tubulin polyglutamylase that forms polyglutamate side chains on tubulin. Probably acts when complexed with other proteins.

Gene ID:

TTLL9

Uniprot

Q3SXZ7

Synonyms:

Probable tubulin polyglutamylase TTLL9 (EC 6. -. -.) (Tubulin--tyrosine ligase-like protein 9), TTLL9, C20orf125


Immunogen:

Recombinant Human Probable tubulin polyglutamylase TTLL9 protein (1-230AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO35982 at dilution of 1:100.