

SLC25A35 Antibody, Biotin conjugated



PACO35993

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SLC25A35

Uniprot

Q3KQZ1

Synonyms:

Solute carrier family 25 member 35, SLC25A35

Immunogen:

Recombinant Human Solute carrier family 25 member 35 protein (120-168AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A