

SP6 Antibody



PACO35994

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Promotes cell proliferation.

Gene ID:

SP6

Uniprot

Q3SY56

Synonyms:

Transcription factor Sp6 (Krueppel-like factor 14), SP6, KLF14


Immunogen:

Recombinant Human Transcription factor Sp6 protein (101-376AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO35994 at dilution of 1:100.