

ANKDD1A Antibody



PACO36002

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ANKDD1A

Uniprot

Q495B1

Synonyms:

Ankyrin repeat and death domain-containing protein 1A, ANKDD1A


Immunogen:

Recombinant Human Ankyrin repeat and death domain-containing protein 1A protein (272-490AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO36002 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO36002 at dilution of 1:100.