

C5orf22 Antibody



PACO36034

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Gene ID:

C5orf22

Uniprot

Q49AR2

Synonyms:

UPF0489 protein C5orf22, C5orf22


Immunogen:

Recombinant Human UPF0489 protein C5orf22 protein (101-420AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: C5orf22 antibody at 8 μ g/ml + Mouse brain tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 50, 21, 19 kDa. Observed band size: 50 kDa.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO36034 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36034 at dilution of 1:100.