

COX19 Antibody



PACO36062

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

May be required for the assembly of mitochondrial cytochrome c oxidase.

Gene ID:

COX19

Uniprot

Q49B96

Synonyms:

Cytochrome c oxidase assembly protein COX19 (hCOX19), COX19


Immunogen:

Recombinant Human Cytochrome c oxidase assembly protein COX19 protein (8-90AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO36062 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).