

CCDC184 Antibody



PACO36086

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

CCDC184

Uniprot

Q52MB2

Synonyms:

Coiled-coil domain-containing protein 184, CCDC184, C12orf68


Immunogen:

Recombinant Human Coiled-coil domain-containing protein 184 protein (1-194AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO36086 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO36086 at dilution of 1:100.