

PACO36130

---

## Product Information

**Size:**

50ug

**Reactivity:**

Staphylococcus aureus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

fabD

**Uniprot**

Q5HGK3

**Synonyms:**

Malonyl CoA-acyl carrier protein transacylase (MCT) (EC 2.3.1.39), fabD

**Immunogen:**

Recombinant Staphylococcus aureus Malonyl CoA-acyl carrier protein transacylase protein (1-308AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A