

# Ribosome-inactivating protein gynostemmin Antibody PACO36174



---

## Product Information

**Size:**

50ug

**Reactivity:**

Gynostemma pentaphyllum

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Ribosome-inactivating protein of type 1, inhibits protein synthesis in animal cells.

**Gene ID:****Uniprot**

Q5F0I3

**Synonyms:**

Ribosome-inactivating protein gynostemmin (EC 3.2.2.22) (rRNA N-glycosidase)

**Immunogen:**

Recombinant Gynostemma pentaphyllum Ribosome-inactivating protein gynostemmin protein (24-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A