

Ribosome-inactivating protein gynostemmin
Antibody, Biotin conjugated
PACO36177



Product Information

Size:

50ug

Reactivity:

Gynostemma pentaphyllum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ribosome-inactivating protein of type 1, inhibits protein synthesis in animal cells.

Gene ID:

Uniprot

Q5F0I3

Synonyms:

Ribosome-inactivating protein gynostemmin (EC 3.2.2.22) (rRNA N-glycosidase)

Immunogen:

Recombinant Gynostemma pentaphyllum Ribosome-inactivating protein gynostemmin protein (24-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A