

IGSF11 Antibody



PACO36198

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Functions as a cell adhesion molecule through homophilic interaction. Stimulates cell growth.

Gene ID:

IGSF11

Uniprot

Q5DX21

Synonyms:

Immunoglobulin superfamily member 11 (IgSF11) (Brain and testis-specific immunoglobulin superfamily protein) (Bt-IGSF) (V-set and immunoglobulin domain-containing protein 3), IGSF11, BTIGSF CXADRL1 VSIG3


Immunogen:

Recombinant Human Immunoglobulin superfamily member 11 protein (23-241AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: IGSF11 antibody at 1.5 μ g/ml. Lane 1: Mouse brain tissue. Lane 2: Jurkat whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 47, 44 kDa. Observed band size: 47 kDa.