

# FAM76B Antibody



PACO36218

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FAM76B

**Uniprot**

Q5HYJ3

**Synonyms:**

Protein FAM76B, FAM76B

**Immunogen:**

Recombinant Human Protein FAM76B protein (97-233AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A