

NAALADL2 Antibody, HRP conjugated



PACO36231

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be catalytically inactive.

Gene ID:

NAALADL2

Uniprot

Q58DX5

Synonyms:

Inactive N-acetylated-alpha-linked acid, c dipeptidase-like protein 2 (NAALADase L2), NAALADL2

Immunogen:

Recombinant Human Inactive N-acetylated-alpha-linked acid, c dipeptidase-like protein 2 protein (143-795AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A