

COBLL1 Antibody



PACO36250

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

COBLL1

Uniprot

Q53SF7

Synonyms:

Cordon-bleu protein-like 1, COBLL1, KIAA0977


Immunogen:

Recombinant Human Cordon-bleu protein-like 1 protein (1-106AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36250 at dilution of 1:100.



Immunofluorescent analysis of Hela cells using PACO36250 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO36250 at dilution of 1:100.