

PDZD11 Antibody



PACO36262

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

PDZD11

Uniprot

Q5EBL8

Synonyms:

PDZ domain-containing protein 11 (ATPase-interacting PDZ protein) (Plasma membrane calcium ATPase-interacting single-PDZ protein) (PMCA-interacting single-PDZ protein), PDZD11, AIPP1 PDZK11 PISP


Immunogen:

Recombinant Human PDZ domain-containing protein 11 protein (1-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human appendix tissue using PACO36262 at dilution of 1:100.