

FNDC3B Antibody



PACO36286

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May be a positive regulator of adipogenesis.

Gene ID:

FNDC3B

Uniprot

Q53EP0

Synonyms:

Fibronectin type III domain-containing protein 3B (Factor for adipocyte differentiation 104) (HCV NS5A-binding protein 37), FNDC3B, FAD104 NS5ABP37


Immunogen:

Recombinant Human Fibronectin type III domain-containing protein 3B protein (5-71AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO36286 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO36286 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).